

Specifications:



INSTANT DENTAL IMAGING

Portable Dental X-ray Solution



Tube Voltage:
3.5mA

Tube Current:
70kV

Exposure Time:
0.01~2.0s

Max. Output Power:
245W

Focal Spot to Skin Distance:
≥200m

Inverter Frequency:
200kHz

Battery Capacity:
300+ Images per Full Charge

Charging Method:
Wireless Charging

Weight:
1.8kg

Compatible Sensors:
iRay, Brigates , Handy



All-in-One Design

Nano Lead Shielding Glass

- Reduces scattered radiation
- Better radiation protection

Low-Dose Portable X-ray

- Lightweight, only 1.9 kg
- 4-layer Lead shielding

Built-in Imaging System

- 4" HD touchscreen
- Capture, view and export images

USB Sensor Port

- Compatible with most intraoral sensors

Wireless Charging

- Charging and storage in one

Intraoral Sensor Holder

Accurate positioning for clear images



Long Battery Life

Ready whenever you need it

60 min Fast Charging

300+ Images per Charge



Easy to Use

No manual settings required



AED Smart Exposure Control

Automatically selects the best exposure settings.

Faster Workflow, Better Patient Experience

Improve Clinical Efficiency



Compact, Portable and Easy

Easy to Share Between Operatories

